

The Labcyte® Echo® Cherry Pick software

The Labcyte Echo Cherry Pick software makes it simple and easy to cherry pick samples, and works seamlessly with all Echo liquid handlers. You choose the volumes, wells, and replicates that best fit your needs. Add standards and blanks to a multi-well plate with either traditional placement or in a randomized fashion to reduce geometry-related analysis errors.

The new software allows you to import pick-lists to find and efficiently transfer compounds that are distributed anywhere among the source plates. It also enables you to consolidate compounds selected from several source plates into a single destination plate.

Fast

- Saves time by efficient transfer
- Executes superior data management
- Maintains survey history

Flexible

- Uses pick lists to define well locations and volumes for transfer
- Customizes reports to match your informatics infrastructure
- Generates output files in XML, CSV, and TXT format
- Transfer from 384- or 1536-well Echo qualified source plates to 96-, 384-, 1536-, or 3456-well assay plates
- Allows the output files to be run specific or destination plate specific

Accurate

- Automates manual steps
- Removes washing steps and cross contaminations

Labcyte Catalog number: Echo-CPick
Requires V2.1 or higher operating software

